

## Portable Medium Exchange System



### Features

- Automates the process of medium exchange at a specific time interval.
- Battery-operated. Usable in an incubator.
- A standard 6-well plate can be used with this system.
- Wells are observable on a microscope.
- Can be remodeled into a perfusion culture system.

### Specifications

Model Number	CEME-0102
Volume of Medium Exchange	Approx. 2.7 ml /cycle
Power Supply	AA (LR6) battery x 2 (for 3 days ~ 7 days)
Wetted Materials	PP, Silicone, Stainless Steel 304
Dimensions	194 x 228 x 140mm

Note: Details including specifications may change without notification.

### TAKASAGO ELECTRIC, INC.

No.5 XingHan Street BLK B# 06-01/08 SuZhou Industric Park SuZhou City JiangSu Province China

Tel +86-512-67610522 Fax +86-512-67610533

E-mail: [info@takasago-elec.net](mailto:info@takasago-elec.net) URL: <http://www.takasago-elec.com.cn/>